

GLUBRAN<sup>®</sup> 2

# Applicators



## Wide range for personalised procedures



The **Drop-by-drop** devices are designed to allow accurate applications of Glubran<sup>®</sup> 2 to take advantage of its adhesive property to replace or reinforce sutures, to cover small areas, for haemostasis or for sealing. The product must be applied in a quantity of about one drop per cm<sup>2</sup> in order to achieve the best result.


The **Dispensing tip** for open surgery features a brush tip for applying a thin layer. Particularly suitable for dispensing Glubran<sup>®</sup> 2 to seal anastomoses and reinforce sutures.

The **Flexible and rigid laparoscopic catheters** are indicated for drip application of Glubran<sup>®</sup> 2 as a haemostatic, adhesive or sealant agent in laparoscopic surgical procedures.

The **nebulizers** with autonomous gaseous propulsion are designed to enhance all the properties of Glubran<sup>®</sup> 2, both in laparotomic and laparoscopic surgical procedures, by forming a thin film of product.

### Nebulizers

- Preparation is quick and easy
- These are independent single-use systems
- They do not require pressure regulators
- They do not connect to operating room air fittings or to CO<sub>2</sub> cylinders
- The inert propellant for nebulisation is non-toxic, biocompatible and has an autonomy of 8 minutes of continuous supply
- The outlet pressure is constant and controlled with no risk of air embolism



**Drop-by-drop control application devices**



G2-DCD-210-8T



G-LLS

**Thin film application devices**



G-DT






**Laparoscopic catheters for drop-by-drop control application**



G2-LPC FLEXIBLE




G2-LPC-RIG RIGID

<b>Nebulizers</b>	<b>APPLICATIONS</b>	<b>DISTANCE FROM TISSUE</b>	<b>NEBULIZING DIAMETER</b>	<b>PROPERTIES</b>	<b>RATIONALE</b>
 <b>NEBULIZER DEVICE</b> G2-NBT		5-7 cm	4 cm	Flexible, for large surfaces, easy pointing, also suitable for open surgery	
 <b>SMALL NEBULIZER DEVICE</b> G2-NBT-SMALL	In laparoscopy thoracoscopy	2-3 cm	1 cm	Flexible, for small surfaces, easy pointing, also suitable for open surgery	SEALING
 <b>RIGID NEBULIZER DEVICE</b> G2-NBT-RIG		5-7 cm	4 cm	Rigid, for large surfaces	HAEMOSTASIS REINFORCEMENT
 <b>SHORT NEBULIZER DEVICE</b> G2-NBT-SHORT		5-7 cm	4 cm	Short, for large surfaces	AEROSTASIS BILIOSTASIS
 <b>SHORT MINI NEBULIZER DEVICE</b> G2-NBT-SM-SHORT	In open surgery	1 cm	1 cm	Short, for small surfaces	

**GLUTÄCK**

**Laparoscopic device for atraumatic fixation of hernia meshes**



GB-DS-25  
GB-DS-50

Available in 2 sizes, 25 or 50 drops, depending on the type of procedure and size of the hernia mesh

