GLUBRAN[®] Tiss 2

CE 0373

COLLE À USAGE CUTANÉ BI-COMPOSANT

NBCA (n- butyle 2 cyanoacrylate) monomère, OCA (2- octyl cyanoacrylate) monomère

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Bibliographie

1. Lombardi CP, Bracaglia R, Revelli L, Insalaco C, Pennestri F, Bellantone R, Raffaelli M. **Aesthetic result of thyroidectomy: evaluation of different kinds of skin suture.** Ann Ital Chir. 2011 Nov-Dec;82(6):449-55
2. Gulalp B, Seyhan T, Gursoy S, Altinors MN. **Emergency wounds treated with cyanoacrylate and long-term results in pediatrics: a series of cases; what are the advantages and boards?** BMC Res Notes. 2009 Jul 14;2:132.
3. Chibbaro S, Tacconi L. **Use of skin glue versus traditional wound closure methods in brain surgery: A prospective, randomized, controlled study.** J Clin Neurosci. 2009 Apr;16(4):535-9.
4. Gennari R, Rotmensz N, Ballardini B, Scevola S, Perego E, Zanini V, Costa A. **A prospective, randomized, controlled clinical trial of tissue adhesive (2-octylcyanoacrylate) versus standard wound closure in breast surgery.** Surgery. 2004 Sep;136(3):593-9.
5. Graefe H, Wollenberg B, Brocks C. **Application of skin adhesives in head and neck surgery: analysis of cosmetic results, applicability and cost-effectiveness of cyanoacrylate-based adhesives.** HNO. 2008 Sep;56(9):947-54.
6. Craven NM, Telfer NR. **An open study of tissue adhesive in full-thickness skin grafting.** J Am Acad Dermatol. 1999 Apr;40(4):607-11.
7. Farion K, Osmond MH, Hartling L, Russell K, Klassen T, Crumley E, Wiebe N. **Tissue adhesives for traumatic lacerations in children and adults.** Cochrane Database Syst Rev. 2002;(3):CD003326. Review.
8. Haghpanah S, Vafafar A, Golzadeh MH, Ardeshtiri R, Karimi M. **Use of Glubran 2 and Glubran tissue skin adhesive in patients with hereditary bleeding disorders undergoing circumcision and dental extraction.** Ann Hematol. 2011 Apr;90(4):463-8
9. Amiel GE, Sukhotnik I, Kowar B, Siplovich L. **Use of N-butyl-2-cyanoacrylate in elective surgical incisions—longterm outcomes.** J Am Coll Surg. 1999 Jul;189(1):21-5.
10. **Topical Skin Adhesive: Strong, Safe, Secure Wound Closure.** Ed. Pfiedler Enterprises 2101 S. Blackhawk St, Suite 220, Aurora, CO 80014 www.pfiedlerenterprises.com
11. Hosseini SM, Rasekhi AR, Zarenezhad M, Hedjazi A. **Cyanoacrylate glue dressing for hypospadias surgery.** N Am J Med Sci. 2012 Jul;4(7):320-2.
12. Howell JM, Bresnahan KA, Stair TO, Dhindsa HS, Edwards BA. **Comparison of effects of suture and cyanoacrylate tissue adhesive on bacterial counts in contaminated lacerations.** Antimicrob Agents Chemother. 1995 Feb;39(2):559-60
13. Amin M, Glynn F, Timon C. **Randomized trial of tissue adhesive vs staples in thyroidectomy integrating patient satisfaction and Manchester score.** Otolaryngol Head Neck Surg. 2009 May;140(5):703-8.
14. Data on GEM S.r.l. file.
15. Singer AJ, Quinn JV, Hollander JE. **The cyanoacrylate topical skin adhesives.** Am J Emerg Med. 2008 May;26(4):490-6.
16. Osmond MH, Klassen TP, Quinn JV. **Economic comparison of a tissue adhesive and suturing in the repair of pediatric facial lacerations.** J Pediatr. 1995 Jun;126(6):892-5.

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TUV America Inc ISO 13485:2003
GEM S.r.l. Via dei Campi 2 - PO Box 427 - 55049 Viareggio (LU) Italy