A.5650

A.5651

A.5660

A.5661

LEM[®]

Legatore Emorroidario Monouso

Disposable Haemorrhoid Banding Ligator

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The LEM Disposable Haemorrhoid Banding Ligator is a unique, patented design made from a light, tough polymer. It is a single-patient product available in two variants: for use with either suction or grasping forceps.

Both variants are supplied with a stepped cylindrical cone (Ø 10mm) for easy band loading, or with a smooth cone and a ring-loader. The 20-degree angled tip of the LEM ensures optimal visual clarity and the ergonomically designed trigger enables smooth and accurate banding procedures. A length of flexible suction tubing connects the LEM to a low-pressure suction aspirator. The haemorrhoid is retracted into the LEM's barrel-end and then banded by the operator using the trigger-action. Throughout the procedure, the suction may be manually thumb-controlled by the operator. If the forceps variant is chosen, suction is replaced by the insertion of grasping forceps through a cutaway section in the barrel so that the haemorrhoid may be retracted manually prior to being banded.

The LEM is designed for use with the Sapi Med Self Light operative proctoscopes.

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
Suction tubing connects the LEM to a low-pressure suction aspirator, and the LEM is then inserted into an operative proctoscope and placed over the surgical site. The haemorrhoid is retracted from the tissue by suction and then banded using the trigger-mechanism. During the procedure, the rate of suction may be controlled manually by the practitioner's thumb. Latex free. Non sterile.

LEM a TRAZIONE - LEM for use with FORCEPS

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In place of suction, grasping forceps are inserted through a cut-out opening in the barrel of the bander and used to manually retract the haemorrhoid for banding. Latex free. Non sterile.

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
REF A.5650: LEM - Disposable haemorrhoid ligation suction and banding instrument.

REF A.5650.3: LEM - Disposable haemorrhoid ligation suction and banding instrument complete with 3 rings.

REF A.5650.5: LEM - Disposable haemorrhoid ligation suction and banding instrument complete with 5 rings.

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
REF A.5651: LEM - Disposable haemorrhoid banding ligator for use with suction, complete with rings loader.

REF A.5651.3: LEM - Disposable haemorrhoid banding ligator for use with suction, complete with rings loader and 3 rings.

REF A.5651.5: LEM - Disposable haemorrhoid banding ligator for use with suction, complete with rings loader and 5 rings.

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
REF A.5660: LEM - Disposable haemorrhoid banding ligator for use with forceps,

REF A.5660.3: LEM - Disposable haemorrhoid banding ligator for use with forceps, complete with 3 rings.

REF A.5660.5: LEM - Disposable haemorrhoid banding ligator for use with forceps, complete with 5 rings.

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REF A.5661: LEM - Disposable haemorrhoid banding ligator for use with forceps, complete with rings loader

REF A.5661.3: LEM - Disposable haemorrhoid banding ligator for use with forceps, complete with rings loader and 3 rings.

REF A.5661.5: LEM - Disposable haemorrhoid banding ligator for use with forceps, complete with rings loader and 5 rings.

Quantity per box: 25 pieces

Distributed by