

Vital Readings within Reach

M3A Vital Signs Monitor



Portable Design



M3A Vital Signs Monitor

12,000

Measurements
Review

EDAN M3A Vital Signs Monitor provides essential measurement of patient's SpO₂, NIBP and Temperature. It offers an effective solution and brings high performance vital signs spot-check at your finger tips.

- 3.5" high resolution color TFT-LCD screen & LED display
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Auto/Manual/Continuous/Average BP mode
- Real-time measurements & trend display
- Built-in rechargeable Lithium-ion battery for 18 hours of continuous work
- Compatible with CMS-Lite for PC data management
- Nurse call
- Suitable for adult, pediatric and neonatal patients
- Barcode scanner support (optional)
- Ethernet connection for bi-directional communication with MFM-CMS central monitoring system (optional)
- Optional Infrared Ear Temp and Quick Temp



Infrared Ear Temp



EDAN T2 Quick Temp



Covidien Filac 3000 Quick Temp

Configuration Table

	SpO ₂	NIBP	Temp
SpO ₂ only	●		
NIBP only		●	
SpO ₂ +NIBP	●	●	
NIBP+ Temp		●	●
SpO ₂ +NIBP+ Temp	●	●	●

Note: The Temp could be any one of RADIANT THP59J, EDAN T2 Quick Temp and Covidien F3000 Quick Temp.

Mounting Solutions



Wall Mount



Rolling Stand



Bed Rail Clamp



A world of potential

Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-PM-M3A
V1.5-20240724