

ART. 678

OSTOMY SUPPORT BINDER

CHARACTERISTICS AND FUNCTIONALITIES

Helps recover as fast and better as possible:

- Adjustable and homogeneous compressive capacity to prevent abdominal hernias

User friendly:

- Replaceable bacteriostatic frontal pannel
- 7 cm diameter decentralized colostomy hole (customizable up to 11 cm) to adjust the binder to any necessity
- Bag cover in elastic and ultra-light cotton
- Unisex
- Completely openable

Comfortable:


- Breathable and multi extensive elastic fabric (Pavis "Airflex" woven)

Hygienic:

- 100% soft cotton in contact with the skin
- Hypoallergenic fabric

INDICATIONS


- Colostomies
- Ileostomies
- Ureterostomies



ART. RCMB678
PANNEL

ANTIBACTERIAL FABRIC

AirGeniX.
ANTIBACTERIAL FABRIC



SIZE	S	M	L
CM	75 - 90	90 - 115	115 - 135

HEIGHT: 22 CM

COLOURS	
---------	--

**100%
COTTON**
ON YOUR SKIN

