


TECHNICAL SPECIFICATIONS



Autonomy

Up to five hours of bright light.

Adaptable

Can be fitted to any surgical instrument with a diameter of 5 to 8mm.

Versatile

Up to 3 RUTILIGHT® devices can be connected to the same power source.

Thrifty

Low power consumption.

Optimization

Optimal colour temperature for organ, vessel and tissue colour rendition.



Port Ginesta – Local 814 – Mod. B
08860 Castelldefels (Barcelona)
Phone: +34 938 325 897
E-mail: info@swanmedical.es
www.swanmedical.es


RUTILIGHT

The light beyond your hands




LIGHTS UP SURGERY

LIGHTS UP SURGERY




The RUTILIGHT® device allows for the comfortable and practical illumination of the surgical field regardless of the surrounding lighting, enabling access to badly lit areas even when using surgical headlights.

The RUTILIGHT® device improves the surgeon's visual field in all circumstances, including deep incisions and those that are difficult to access, minimising complication rates, facilitating surgery and reducing operation time.



MAXIMUM VERSATILITY



The RUTILIGHT® device is a light ring featuring LED technology which can be adapted to any surgical instrument, in the electroscalpel handle or in any other necessary one in a surgery, without affecting function.

FOR MOST SURGERIES

Coloproctology
Colon, rectal, anal and endoanal surgery.

Endocrinology
Thyroid, parathyroid and suprarenal surgery and ganglion drainage.

Breast surgery
Subcutaneous areola-sparing mastectomy and underarm surgery.

Oesophagogastric surgery
In the abdominal and cervical phases.

General and digestive surgery
Hepato-bilio-pancreatic and other surgeries.

Pediatric surgery
Solid tumours, severe malformations and stenosis.

Plastic and reconstructive surgery
Complex maxillofacial and other operations.

Morbid obesity
In all kind of open surgeries.

Transplants
Liver, kidney, pancreas and lung.

Peritoneal oncology
Peritoneal carcinomatosis.

Traumatology
Complex shoulder, hip, knee, foot and ankle injuries.

Urology
Pathologies affecting the urinary tract and reproductive system.

Gynecological
Hysterectomy, abdominal myomectomy and complex surgeries.

REFERENCES

RUTILIGHT REF.: RL-01

Each box contains 25 RUTILIGHT® devices in sterile single-use individual packages.

POWER SOURCE REF.: FA-RL01

