

Reusable

Soft and comfortable material

Latex free



Bi-Mask

Anatomically shaped anaesthesia face mask

/ Bi-Mask

The Bi-Mask is a face mask made of silicone with integrated plastic shell. This supports a comfortable grip, as the mask is malleable and yet provides stability without collapse. Due to the anatomical shape, the Bi-Mask provides an airtight seal with minimum pressure and the ergonomic design allows an optimised grip via C-grip. The Bi-Mask is available in 5 sizes and processable up to 100 times (autoclavable at 134°C).

Features



Order information

Bi-Mask / Reusable

Size	Patient	Connection	Colour	REF	Box
2	Child, small	22 mm I.D.	■	35-65-222	1
3	Child, large	22 mm I.D.	■	35-65-223	1
4	Adult, small	22 mm I.D.	■	35-65-224	1
5	Adult, medium	22 mm I.D.	■	35-65-225	1
6	Adult, large	22 mm I.D.	■	35-65-226	1



Hook Ring / Reusable

REF	Box
88-60-266	10

Fixation Tape / For single use

Patient	REF	Box
Adult	35-30-255	10

Fixation Tape / Reusable

Patient	REF	Box
Child	35-70-155	10
Adult	35-70-255	10

This product is manufactured without the use of natural rubber latex, unless otherwise specified.
This product does not contain phthalates which require labelling according to CLP Regulation (EC) 1272/2008.